VVP

DEUTSCH

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ENGLISH

VVP is a solution of 7% or 10% pharma grade Polyvinylpyrrolidone (PVP) in HEPES buffered Human Tubal Fluid (HTF) Medium with Human Serum Albumin (HSA), ready to use. VVP is used for reducing the motility and to simplify the immobilization and easy handling of sperm during intracellular microinjection (ICSI). VVP also prevents sperm sticking at the injection pipette (VICISIPET) and gives better control of media and sperm in the ICSI pipette.

COMPOSITION

- Polyvinylpyrrolidone
- Human Serum Albumin
- Physiological Salts
- Pyruvate
- Lactate
- Glucose
- HEPES HTF

PRODUCT SPECIFICATIONS

- pH: between 7.25 - 7.45
- Sterility: sterile [SAL 10-3]
- Mouse Embryo Assay (Blastocysts after 9h): ≥ 80%
- A certificate of analysis and MSDS are available upon request.
- Endotoxins: < 0.5 EU/ml

STERILE A Sterilized using aseptic processing techniques.

PRE-USE CHECK

- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if any discoloration, turbidity or precipitation is noticed.
- Do not use the product if seal of the container is defect or opened upon delivery.

INSTRUCTIONS FOR USE

- ① Place several droplets of HEPES buffered media (5 µl to 10 µl) and a droplet of VVP (5 µl to 10 µl) in a petri dish, cover with VOLL and warm up to 37°C.
- ② Add 1 µl to 2 µl of washed sperm to the VVP droplet.
- ③ Immobilize and aspirate the sperm for injection into the oocytes during ICSI.

STORAGE AND STABILITY

- Store VVP between 2°C to 8°C.
- Do not freeze before use.
- Leave product in original packaging and protect from light.
- The product is stable until the expiration date shown on the label or within 7 days of the date of first opening, provided that the proper aseptic procedures have been achieved by the user.

WARNING AND PRECAUTIONS

VITROMED recommends handling of media only with aseptic techniques. Best measures to prevent infections resulting from use of medical products prepared from human blood or plasma include precise screening of each donation and plasmapools for distinct markers of infection and the use of validated methods for inactivation of viruses during production. Despite this, a possible transmission cannot completely ruled out. To the current state, there are no reports about confirmed cases about transmissions of viral infections or Creutzfeldt-Jakob disease (CJD) due to human albumin. Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (LAF-bench ISO class 5) to avoid contamination. Ensure a correct disposal.

DE / EN

Artikelbezeichnung Name of product	VVP
REF	V-PVP-7 V-PVP-10
Stück pro VE	PVP 7% 5 x 0,5 ml Phiole PVP 10% 5 x 0,5 ml Phiole
Piece per unit	PVP 7% 5 x 0,5 ml vial PVP 10% 5 x 0,5 ml vial

Technischer Support
For technical support

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VITROMED GmbH