## DEUTSCH

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## ENGLISH

SPERM FREEZE is a ready-to-use hepes buffered cryopreservation medium for In Vitro Fertilization (IVF) and Intracellular microinjection (ICSI) and used for cryopreservation and storage of human sperm.

### COMPOSITION

- HEPES
- Lactate
- Sucrose
- Glycine
- Glycerol
- Phenol Red
- Human Serum Albumin (4 g/L)
- Physiological Salts
- Dextrose Monohydrate

### PRODUCT SPECIFICATIONS

- pH: between 7.20 - 7.90 (Release criteria: 7.20 - 7.60)
- Endotoxins: < 0.25 EU/ml
- Sterility: sterile (SAL 10<sup>-3</sup>)
- Human Sperm Survival Assay tested ≥ 80% survival after 4 h exposure of untreated semen to the test medium
- A certificate of analysis and MSDS are available upon request.

**STERILE A** Sterilized using aseptic processing techniques.

### PRE-USE CHECK

- Do not use the product if any signs of microbial contamination are noticed.
- Do not use if any discoloration, turbidity or precipitation is noticed.
- Do not use the product if seal of the container is defect or opened upon delivery.

### INSTRUCTIONS FOR USE

- Before use warm up to room temperature and mix it thoroughly.
- Remove the desired volume of product using aseptic techniques.
- Write the date of first opening on the product label.
- Once removed, do not return any volume of product to the original container.
- SPERM FREEZE should only be used for its intended purpose by qualified and trained personnel only.

### FREEZING

1. Prior to the freezing procedure allow the semen to liquify for 30 minutes at room temperature.
2. Mix 1 ml of sperm with 0.7 ml of SPERM FREEZE. The medium should be added drop by drop into the semen and the suspension carefully mixed after each addition. Caution: To avoid cold-shock, make sure SPERM FREEZE is at room temperature.
3. Leave the mixture at room temperature for 10 minutes.
4. Load the cryo vials or straws with the sperm suspension according to the instructions of the manufacturer.
5. For straws mind the position of air bubbles according to the instructions of the manufacturer.
6. Place the vials or straws for 15 minutes vertically in a Dewar vessel just above the level of the liquid nitrogen.
7. Plunge the straws or cryo vials completely into the liquid nitrogen and transfer them to a storage tank.

### THAWING

1. After removing the required straws or cryo tubes from the liquid nitrogen, warm the straws at room temperature for 5 minutes.
2. Open the straws or cryo tubes following the manufacturer's instructions and transfer the thawed sperm suspension in a test tube.
3. Dilute the sperm suspension with SPERM WASH (3 mL per 0.5 mL semen), mix properly and centrifuge 15 minutes at 300 x g to 350 x g.
4. Resuspend the pellet in SPERM WASH.
5. Concentrate the sperms by a density gradient centrifugation or swim-up procedure by the use of GRAD and SPERM WASH.

### STORAGE AND STABILITY

- Store SPERM FREEZE between 2°C to 8°C.
- Do not freeze before use.
- Keep away from sunlight.
- The product is stable until the expiration date shown on the label or within 7 days of the date of first opening, provided that the proper aseptic procedures have been achieved by the user.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).

### WARNING AND PRECAUTIONS

VITROMED recommends handling of media only with aseptic techniques. Best measures to prevent infections resulting from use of medical products prepared from human blood or plasma include pre-screening of each donation and plasmapools for distinct markers of infection and the use of validated methods for inactivation of viruses during production. Despite this, a possible transmission cannot completely ruled out. To the current state, there are no reports about confirmed cases about transmissions of viral infections or Creutzfeldt-Jakob disease (CJD) due to human albumin. Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (LAF-bench ISO class 5) to avoid contamination. Ensure a correct disposal.

## DE / EN

Artikelbezeichnung Name of product	<b>SPERM FREEZE</b>
REF	V-SFM-12
Stück pro VE Piece per unit	4 x 12 ml Glas-Flasche 4 x 12 ml glass bottle

Technischer Support  
For technical support

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**VITROMED GmbH**